PATHOZYME® PANCREATIC/GUT CANCER ANTIGEN 19-9 Ref OD277 Enzimski imunotest za kvantitativno određivanje gastrointestinalnog CA19-9 u humanom serumu Čuvati na temperaturi od 2-8 °C. NE ZAMRZAVATI. Samo za in vitro upotrebu.

OPREZ

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Testovi daju dobru korelaciju.

REFERENCE

1. Glenn, J., Steinberg, W.M., Kurtzman, S.H., et al. Evaluation of the utility of a radioimmunoassay for serum Cancer Antigen 19-9 level in patients before and after treatment of carcinoma of the pancreas. J. Clin. Oncol. 1988;6:462-8.
2. Hayakawa, T., Kondo, T., Shibata, T. et al. Sensitive serum markers for detecting pancreatic cancer. Cancer 1988;61:1827-31.
3. Koprowski, H., Herly, Steplewski, Z., et al. Specific antigen in serum of patients with colon carcinoma. Science 1981;212:53-55.
4. Malesci, A., Tommasini, M.A., Bonato, C. et al. Determination of Cancer Antigen 19-9 antigen in serum and pancreatic juice for differential diagnosis of pancreatic adenocarcinoma from chronic pancreatitis. Gastroenterology 1987;92:60-7.
5. Safi, F., Roscher, R., Bittner, R., et al. High sensitivity and specificity of CA 19-9 for pancreatic carcinoma in comparison to chronic pancreatitis. Serological and immunohistochemical findings. Pancreas 1987;2:398-403.
6. Steinberg, W. The clinical utility of CA 19-9 tumour associated antigen. American J. of Gastroenterology 1990;85:350-355.
7. Steinberg, W.M., Gelfand, R., Anderson, K.K., et al. Comparison of the sensitivity and specificity of the CA 19-9 and Carcinoembryonic antigen assays in detecting cancer of the pancreas. Gastroenterology 1986;90:343-9.
8. Takasaki, H., Uchida, E., Tempero, M.A., et al. Correlative study on expression of CA 19-9 and DU-Pan-2 in tumour tissue and in serum of pancreatic cancer patients. Cancer Res. 1988;48:1435-8.
9. Tetsuta, M., Yamamura, H., Iishi H., et al. Values of Cancer Antigen 19-9 in the serum, pure pancreatic juice and aspirated pancreatic material in the diagnosis of malignant pancreatic tumour. Cancer 1985;56:2669-73.
10. Wang, T.H., In Chen, J.W., et al. Non-invasive diagnosis of advanced pancreatic cancer by real-time ultra sonography. Carcinoembryonic antigen and carbohydrate antigen 19-9. Pancreas 1986;1:1219-23.
11. Strom B.L., Maislin, G., West, S. I. et al. Serum CEA and Cancer Antigen 19-9; potential future diagnostic or screening tests for gallbladder cancer? Int. J. Cancer 1990;45:821.

BRZI POSTUPAK ISPITIVANJA

1. Preneti 10 μ l uzorka ili standarda i 100 μ l pufera za ispitivanje u svaku čašicu i blago mućkati 30 sekundi.
2. Inkubirati 90 minuta na 37 °C.
3. Odbaciti sadržaj čašica i ispirati 5 puta destilovanom vodom.
4. Dodati 100 μ l radnog rastvora konjugata u svaku čašicu. Blago mućkati 10 sekundi.
5. Inkubirati 90 minuta na 37 °C.
6. Odbaciti sadržaj čašica i ispirati 5 puta destilovanom vodom.
7. Dodati 100 μ l rastvora supstrata u svaku čašicu. Blago mućkati 10 sekundi.
8. Inkubirati na tamnom mestu 20 minuta na sobnoj temperaturi (20 °C to 25°C).
9. Dodati 100 μ l rastvora za zaustavljanje reakcije i blago mućkati 30 sekundi.
10. Odmah očitati optičku gustinu (ne nakon isteka 10 minuta) na 450nm.

8077 ISSUE 8 Revised February 2010. **SERBIAN**
©Omega Diagnostics Ltd 2010.

OMEGA DIAGNOSTICS LTD.
Omega House, Hillfoots Business Village
Alva FK12 5DQ, Scotland, United Kingdom
odl@omegadiagnostics.co.uk
www.omegadiagnostics.com
AN ISO 9001 and ISO 13485 CERTIFIED COMPANY